

# Accessories

RADIO NURSE CALL SPECIALISTS

**SLS** | Secure Location  
Solutions Ltd.

medicare  
SYSTEMS



# Assistive Technology

## **MED 30 Rechargeable EPOC Pager and Holster**

Robust, rechargeable pager to suit any care or hospital environment. Complete with docking station for recharging. Receives calls from the nurse call system via a paging transmitter. Pagers are used to alert staff that are anywhere on site. Has a colour coded LCD display screen for easy recognition of call priority.

Power requirements.

Normal operating conditions, battery life 8-12 hours.

Docking station – standard socket outlet 240v a.c.



## **MED 22 PIR Movement Detector**

The PIR movement sensor detects movement in the room. The device connects to the Medicare call point via socket 1. Ideal for residents who are unsteady on their feet.

The PIR sensor is battery operated, a power supply is available (MED 23 PIR Movement Detector PSU).

Complete with a HTM 6022 Accessory Call Adapter.



## **MED 24 Bed Weight Sensor**

The sensor is triggered when the user leaves the bed, releasing the pressure on the unit and activating an accessory call. The weight sensor is used with the Medicare call point using socket 1.

Complete with a HTM 6021 Bed/Chair Adapter. (Max weight 200 kg).



# Assistive Technology

## **MED 11 Epilepsy Detection Unit**

Monitors bed movement, sound, incontinence, vomiting and bed vacation. The movement detection has the capability of supporting those with tonic only (rigid) seizures.

Also available, MED 12 Epilepsy Detection Unit with Cotton Enuresis Sheet, which will trigger an accessory call when the mat comes into contact with moisture. The cotton sheet can be washed at normal temperatures once the connecting cables are disconnected.



---

## **MED 13 Enuresis/Incontinence Detector Kit – Controller with Sheet**

The sheet goes over the mattress and plugs into the room unit via the detection controller and P145 & P141E leads. The call unit will trigger an accessory call when the mat comes into contact with moisture. The cotton sheet can be washed at normal temperatures once the connecting cables are disconnected.

Lasts up to 60 washes.



---

## **MED 14 Washable Cotton Enuresis Sheet Standard**

Replacement enuresis detection sheets for use with MED 12 and MED 13. Standard size, lasts up to 60 washes.

## **MED 15 Washable Enuresis Sheet Large**

Replacement enuresis detection sheets for use with MED 12 and MED 13. Large size, lasts up to 90 washes.



# Pressure Mats

## **MED 7 Standard Economy Floor Pressure Mat**

The economy floor mat can be used near a bed to alert if a resident has left their bed. It can also be used near a door to alert if a resident is about to leave their room.

When pressure is applied to the mat, it triggers an accessory call on the display panel. This should be placed under a carpet or mat.

Complete with a HTM 6022 Floor Mat Adapter.



## **MED 8 Heavy Duty Pressure Mat**

The MED 8 is a heavy duty, carpeted, non-slip mat for use near the resident's bed or door. The mat triggers an accessory call when pressure is applied.

Complete with a HTM 6022 Floor Mat Adapter.



## **MED 9 Chair Pressure Mat**

For use with a HTM6000 or HTM6500 series call point which automatically monitors pressure applied to the chair mat and activates an accessory call when the pressure is released.

Complete with a HTM 6021 Bed/Chair Adapter.



## **MED 10 Bed Pressure Mat**

For use with a HTM6000 or HTM6500 series call point which automatically monitors pressure applied to the bed mat and activates an accessory call when the pressure is released.

Complete with a HTM 6021 Bed/Chair Adapter.





# Assistive Switches

## MED 26 Jelly Bean Switch

An extremely sensitive switch which activates when tapped. Ideal for those with limited strength, movement or eye sight. This switch is available in a range of different colours with interchangeable tops, and the option of two different sizes.

Complete with a HTM 6024 Normal Call Adapter.



## MED 27 Finger Switch

Durable finger mounted switch, ideal for individuals with limited movement. The switch simply wraps around the thumb or finger with a button inside, this can be pressed using the thumb or by pressing the wrapped finger onto a surface. Once pressed, the switch will activate a call on the display panel.

Complete with a HTM 6024 Normal Call Adapter.



## HTM 6025 Blow Switch – with Adapter

This specialist switch is ideal for those who have very limited movement and allows the user to suck or blow into the mouth piece to activate a call on the display panel.

Complete with a HTM 6024 Normal Call Adapter.



## MED 41 Head Tilt Switch

A switch worn by the user as a band around the head and operated by tilting forwards or sideways. Ideal for those with limited movement.

Complete with a HTM 6024 Normal Call Adapter.



# The Choice of The Professionals

## **Assistive Technology tailored to meet residents' individual needs**

At Medicare Systems we understand the importance of giving residents the ability to use a nurse call system, however for some pressing a button is simply not possible.

The Medicare in-house development team explores any new advances in assistive technology that are likely to have a positive involvement in Care Homes and Hospitals.

We manufacture and supply a range of equipment. The Medicare team have come up with innovative solutions for residents with the most challenging conditions.

Medicare Systems look at the needs of individual residents and are happy to discuss any requirements.



6-9 Trinity Street, Dublin 2.

**Phone:** +353 1 617 7887

**Email:** [info@slslimited.ie](mailto:info@slslimited.ie)

**Web:** [www.slslimited.ie](http://www.slslimited.ie)



medicare  
SYSTEMS

[www.medicareystems.co.uk](http://www.medicareystems.co.uk)